


Effector Array: data and analysis

Margie Ackerman
2/21/13

Luminex Bead Assay



Current Status (Antigens)


HIV gp120s	Other Env	Other HIV	Flu Antigens	Other Viral	Controls
gp120 YU2 (37) gp120 Chiang Mai (14) gp120 from Ecoli(54) gp120 93TH975 (19) gp120 BAL (65) gp120 IIIB (CHO) (25) gp120 PVO (35) gp120 CM235 (77) gp120 CM244 (78) gp120 JRCSF (63) gp120 MN (72) gp120 Du151 (73) gp120 SF162 (74) gp120 ZM109F (76) gp120 TRO (61) gp120 Du156.12 (62) gp120 RSC3 (55)	gp41 (HXBc2) (27) gp140 (HXBc2) (52) gp140 CN54 (39) gp140 Du151 (75) SOSIP (45)	p24 (HXBc2) (26) p24 IIIB (33) HIV1 Integrase (28) HIV1 Nef (29) HIV1 Rev (wt) (36) Vif (Ecoli) (38) HIV pr55 Gag (20)	Flu MNA (13) Flu N2 (15) Flu BNA (18) HA Brisbane/10/07 (44) HA NewCal/20/99 (42) HA Wisconsin/67/05 (43) HA (Δ TM) Perth/16/09 (1) HA (Δ TM) Florida/06 (2) HA1 SolomonIsl/06 (3)	HAV P2C-P3A (4) HSV-1 gG (5) HCV E2 (6) PV1 (polio vaccine) PV2 PV3	ahulG (34) ahulGM (12) ahulGA (51) tetanus Ag (54)

~50 Ag; currently working on whole HIV virus bead sets


Current Status (detection)

Subclassing Reagents	FcγRs	Complement Proteins	Lectins
anti-IgG1	FcγRIIa	C1q	SNA
anti-IgG2	FcγRIIb	MBL	WGA
anti-IgG3	FcγRIIIa	C3b	ConA
anti-IgG4	FcγRIIIb	C4b	GNL
	FcγRI		MAL1
	FcRn		LCA
			RCA
Key			PNL
Works on all bead sets			AAL
Works on some bead sets			VVL
Some issues			UAE1
Currently not working			






Validating FcgR Data



- Fc-mutant B12 antibodies from Brian Moldt (Burton Lab - Scripps).
- Compare published ELISA binding to Luminex MFI.




Clinical samples exhibit differential FcγR binding







Brown et al, in preparation

Validating lectin binding




Lectin signals on array correlate with HPLC glycan data


Lectin binding differs among subjects





Similarities among features



Similarities among subjects




Subject-subject correlations, clustered




correlation coefficient

Their antigen-specific MFIs



Feature : function relationships




Features seem predictive

Antibody features

features

subjects




antigens

- gp120
- gp140
- gp41
- HIV1 IIIb Gag
- HIV1 Integrase
- HIV1 Vif (Ecoli)
- p24

readouts





- C1q
- FcγRIIa
- FcγRIIb
- FcγRIIIa
- IgG1
- IgG2
- IgG3
- IgG4

normalized MFI




Clusters driven by
combination of read-
out and antigen

Another view



Clusters driven by
combination of read-
out and antigen

Particular feature relationships (for a particular gp120 antigen)




Correlations among FcγR, IgG1&3; bimodality


Supervised learning

training:


features



& corresponding functions



new subject



(subset of features)




predicted function




Lasso / logistic regression / random forest / support vector machine / Gaussian process

Classification: high vs. low function




Classification example results

Regularized random forest; similar results with logistic- and kernel-based methods



Models effectively use Ab features to qualitatively characterize functions

Regression: quantitative function




repeat for multiple different splits of training/data data

Regression example results


Lasso; similar results with tree- and kernel-based methods (LOOCV)

RFADCC




$r = 0.37$

ADCP



$r = 0.65$

CD107a.MIP1B.IFN γ



$r = 0.68$

Models effectively use Ab features to quantitatively characterize functions